

Transmission Electron Microscope (TEM) sample preparation

There is no perfect method for TEM preparation due to combinations of factors (fixative, buffer, concentration, temperature and processing time). A basic procedure to start TEM work is shown below. However, this should be regarded as a starting point, and methodologies may vary by sample in order to achieve better results. The most important thing to remember is the health of cell(s). It is my experience even that conditions such as time of day and position in cell cycle are important when considering TEM preparation when using synchronized cultures.

The followings are the factors one can change:

Buffer: sodium cacodylate buffer, phosphate buffer, collidine buffer

Fixation method: double fixation, simultaneous fixation

Concentration of fixative: 0.5-5% for glutaraldehyde, 0.1-2% for osmium tetroxide

Duration for fixing: 30 min to 24 hours for glutaraldehyde and osmium tetroxide


Osmic control: with or without sucrose or sorbitol

Concentration of osmotic control: 0.1-0.2M for sucrose and 0.2-0.4M for sorbitol.

Temperature: on ice or room temperature for fixing process

Method: Glutaraldehyde and Osmium Tetroxide (OsO₄) double fixation

- Cool down cell suspension and a centrifuge beforehand.
- Prepare the fixative solution right before processing.
- Use a glass or polypropylene container. Solvents (acetone and propylene oxide) dissolve polystyrene products.



Dysnectes brevis




Rictus lutensis



Wobblia lunata



Caecitellus parvulus




Calkinsia aureus



Cafeteria roenbergensis

Transmission Electron Microscope (TEM) sample preparation



Spurr's resin (Spurr 1969 J. Ultrastr. Res. 26:31) is made from four kinds of monomers and has to be mixed completely. The unpolymerized resin can be kept in a freezer. It needs to be stored separately from water. When removing resin from freezer storage, keep bottle at room temperature without opening the lid. Otherwise condensation forms on the bottle interior. This could be detrimental to the resin.

