

## High salinity medium (for 150; make stock at 300ppt)

(Park 2012 *Extremophiles*)

Medium V (g/L)	150ppt	300ppt
NaCl	137.6	275.2
KCl	3.8	7.6
MgCl <sub>2</sub> •6H <sub>2</sub> O	13.45	26.9
MgSO <sub>4</sub> •7H <sub>2</sub> O	1.65	3.3
CaCl <sub>2</sub> •2H <sub>2</sub> O	0.65	1.3

Fill with RO water/dH<sub>2</sub>O to 1L line, **not** add 1L (~900mL instead)

Uses:

100G and 150G (G = autoclaved barley grain)

For *Trimyema* (microaerophile) -- 12mL 150ppt + 2 barley grains

## ATCC medium: 802 Sonneborn's *Paramecium* medium

*Rye grass Cerophyll*:

Cerophyll\* .....2.5 g

Distilled water.....1.0 L

Add cerophyll to distilled water and boil for 5 minutes. Add 100 ml distilled water to compensate for evaporation. Filter through Whatman #1 filter paper and add 0.5 g Na<sub>2</sub>HPO<sub>4</sub>. Autoclave for 15 minutes at 121C.

SW + 1%LB (SL)

SW + 3%LB (24mL LB in 800mL autoclaved SW)



## Medium TZ

(based on [Mesbah et al. 2007 IJSEM](#))

TZ160 (soda lake medium; quite caustic, wear gloves!) [based on Mesbah et al. 2007 IJSEM]

in g/L; fill to 1L (not “add 1L!”) -- make 1L mark on flask first

KH <sub>2</sub> PO <sub>4</sub>	0.2
MgCl <sub>2</sub>	0.1
KCl	0.2
NH <sub>4</sub> Cl	0.5
NaCl	100
Na <sub>2</sub> CO <sub>3</sub>	68
NaHCO <sub>3</sub>	38

Comes out to ~160-170ppt, pH 9-10, dilute to 20 ppt (measuring with refractometer)

Pour 10mL into 50mL tissue culture flask on its side with a handful (~5-8ish) autoclaved wheat grains as substrate for bacteria

## Resazurin

Stock: 0.2% w/v in dH<sub>2</sub>O

final concentration to be 0.0002% w/v, stock is thus 1000x

filter sterilised and aliquoted into eppies, kept in fridge, in foil (dark)

<https://link.springer.com/article/10.1023/A:1023674315047>

## F/2 Medium

### Directions

For 1 L Total

1. To approximately 950 mL of non-pasteurized seawater ( 30-35ppt), add each of the components in the order specified (except vitamins) while stirring continuously.

2. Bring total volume to 1 L with non-pasturized seawater.

\*For 1.5% agar medium add 15 g of agar into the flask; do not mix.

3. Cover and autoclave medium.

4. When cooled add sterile vitamins.

\*For agar medium add vitamins, mix, and dispense before agar solidifies.

5. Store at refrigerator temperature.

#	Component	Amount	Stock Solution Concentration	Final Concentration
1	NaNO <sub>3</sub> (Fisher BP360-500)	1 mL	7.5 g/100 mL dH2O	880 µM
2	NaH <sub>2</sub> PO <sub>4</sub> ·H <sub>2</sub> O(MCIB 742)	1 mL	0.5 g/100 mL dH2O	36 µM
3	Na <sub>2</sub> SiO <sub>3</sub> ·9H <sub>2</sub> O (Sigma 307815)	1 mL	3 g/100 mL dH2O	106 µM
4	Trace Metals Solution	1 mL/L		
5	Vitamin B <sub>12</sub>	1 mL/L		
6	Biotin Vitamin Solution	1 mL/L		
7	Thiamine Vitamin Solution	1 mL/L		

### Trace Metals Solution

## Directions

For 1 L Total

1. Begin stirring and heating approximately 950 mL of dH<sub>2</sub>O.
2. When the stirring water approaches the boiling point, add the following components in the order listed.
3. Bring the volume to 1 L with dH<sub>2</sub>O.
4. Store at refrigerator temperature.

#	Component	Amount	Stock Solution Concentration	Final Concentration
1	ZnSO <sub>4</sub> ·7H <sub>2</sub> O (Sigma Z 0251)	23 mg/L		0.08 µM
2	MnSO <sub>4</sub> ·H <sub>2</sub> O (Sigma M8179)	152 mg/L		0.9 µM
3	Na <sub>2</sub> MoO <sub>4</sub> ·2H <sub>2</sub> O (J.T. Baker 3764)	7.3 mg/L		0.03 µM
4	CoSO <sub>4</sub> ·7H <sub>2</sub> O (Baker 1696)	14 mg/L		0.05 µM
5	CuCl <sub>2</sub> ·2H <sub>2</sub> O (Fisher C-455)	6.8 mg/L		0.04 µM
6	Fe(NH <sub>4</sub> ) <sub>2</sub> (SO <sub>4</sub> ) <sub>2</sub> ·6H <sub>2</sub> O (Sigma F-1513)	4.6 g/L		11.7 µM
7	Na <sub>2</sub> EDTA·2H <sub>2</sub> O (Sigma ED255)	4.4 g/L		11.7 µM

## WC Medium

## Directions

For 1 L Total:

1. To approximately 900 mL of dd-H<sub>2</sub>O, add the following components in the order listed (not including the vitamins) while stirring continuously.
2. Adjust the pH to 7.8.
3. Transfer the contents of the beaker to a 1-Liter graduated cylinder and bring the total volume to 1.0 Liter with dd-H<sub>2</sub>O.
4. Add the sterile vitamin components and mix well.
5. **Do not autoclave medium! Filter sterilize only.**
6. Store at refrigerator temperature.

#	Component	Amount	Stock Solution Concentration	Final Concentration
1	NaNO <sub>3</sub> (Fisher BP360-500)	1 mL/L	85.1 g/L	1 mM
2	CaCl <sub>2</sub> ·2H <sub>2</sub> O (Sigma C-3881)	1 mL/L	36.76 g/L	0.25 mM
3	MgSO <sub>4</sub> ·7H <sub>2</sub> O (Sigma 230391)	1 mL/L	36.97 g/L	0.15 mM
4	NaHCO <sub>3</sub> (Fisher S 233)	1 mL/L	12.6 g/L	0.15 mM
5	Na <sub>2</sub> SiO <sub>3</sub> ·9H <sub>2</sub> O (Sigma 307815)	1 mL/L	28.42 g/L	0.1 mM
6	K <sub>2</sub> HPO <sub>4</sub> (Sigma P 3786)	1 mL/L	8.71 g/L	0.05 mM
7	H <sub>3</sub> BO <sub>3</sub> (Baker 0084)	1 mL/L	24 g/L	0.39 mM
8	WC Trace Elements Solution	1 mL/L		
9	Vitamin B <sub>12</sub>	1 mL/L		
10	Thiamine Vitamin Solution	1 mL/L		
11	Biotin Vitamin Solution	1 mL/L		

From Lise Forget in Franz Lang's group: (sent 29.09.14)

## WCG and WCL medium, for many protists

WCL medium seems to give better growth than WCG medium, in most cases.

To one liter of distilled and autoclaved water add:

	WCG	WCL	
K <sub>2</sub> HPO <sub>4</sub> *	( 8.71 g/l)	0.25 ml	----
CaCl <sub>2</sub> x2H <sub>2</sub> O	(36.76 g/l)	0.1 ml	1.0 ml
MgSO <sub>4</sub> x7H <sub>2</sub> O	(36.97 g/l)	0.25 ml	1.0 ml
NaNO <sub>3</sub> *	(75.0 g/l)	0.29 ml	1.0 ml
NaHCO <sub>3</sub> *	(12.6 g/l)	0.05 ml	1.0 ml
H <sub>3</sub> BO <sub>3</sub> *	( 6.0 g/l)	0.25 ml	0.5 ml
KCL *	( 7.45 g/l)	----	1.0 ml
NaH <sub>2</sub> PO <sub>4</sub> xH <sub>2</sub> O	( 5.0 g/l)	----	1.0 ml
NH <sub>4</sub> CL *	( 2.65 g/l)	----	1.0 ml
Na <sub>2</sub> EDTA*	(100mM Stock)	0.05 ml	----
Fe/EDTA/trace metals from F/2 stock*		0.25 ml	1.0 ml
F/2 vitamins*		0.5 ml	1.0 ml
(see on F/2 page)			

In order to obtain high cell yields, buffering with HEPES is strongly recommended,

**HEPES buffer:** stock of 1M, adjusted to pH 7.5, end concentration in medium should be about 10-20 mM.

----

### # Misc notes

FeCl<sub>3</sub> (ferric citrate, oxidised state) stock: ~400mM 10<sup>4</sup>x: 3g/50mL dH<sub>2</sub>O, filter sterilised to 1000x in sterile water and then diluted to 100x working solution (appears to be shelf-stable, the higher concentration stocks precipitate), add 200uL / 12mL of culture